



DIZG

DEUTSCHES INSTITUT FÜR
ZELL- UND GEWEBEERSATZ

Gemeinnützige Gesellschaft mbH

DBM pastös

The human demineralised
bone matrix in putty form



DBM pastös – the human demineralised bone matrix in putty form


DBM pastös is an allograft for the treatment of bone defects. The putty is fully converted into new bone after approximately six to twelve months following treatment.

WHAT DOES DBM PASTÖS CONSIST OF?

DBM pastös consists of demineralised bone matrix (DBM) gained from human donors, and a bio-compatible carrier substance. During the production process, the mineral components are first removed from cortical bone. The thus gained demineralised bone matrix is then mixed with a biocompatible substance, a sodium hyaluronate solution.

Hyaluronic acid is present in human tissue such as skin, muscles and bone and lends the bone matrix a user-friendly putty form.

Since 2007, DIZG has provided over 35,000 units to physicians and so helped tens of thousands of patients to improved bone healing.



DBM pastös is approved as a medicinal product in Germany*

SCOPE OF APPLICATIONS

Orthopaedics/traumatology

- > Treatment of fractures in all extremities
- > Filling of defects (ideal for the treatment of solitary bone cysts)
- > Pseudarthrosis

Spinal column

- > Vertebral fusion
- > Cage filling
- > Treatment of fractures

Oromaxillofacial surgery


- > Filling of cysts

Fig 1:

Treatment of bone cysts in OMF surgery



Preoperative X-ray of a follicular cyst



6 months postoperative X-ray



3 years postoperative X-ray

A study from the Charité shows: DBM pastös significantly improves bone regeneration after the removal of cysts in the maxillary region. A total of 90 patients were included in the study, with the cysts showing an average area of 11.2 cm² on the X-ray. The cysts were removed by surgery in all study subjects. In 50 patients, the resulting bone cavities were filled with

human demineralised bone matrix, 40 patients did not receive a filling. The radiological follow-up for a period of up to ten years showed clearly faster new bone formation in the DBM pastös group. The scientists at the Charité explained this effect with the excellent osteoneogenic properties due to the demineralisation of the bone matrix of DBM pastös.⁹

*Marketing authorisation number: PEI.H.03358.01.1

Properties

- > Human demineralised bone matrix
- > Volume ratio of demineralised bone matrix at approx. 76% (v/v)
- > Cell-free⁵
- > Sterile, safe
- > Putty, excellent mouldability (offers perfect sealing with defect margins, can be applied to every form of defect)
- > Bonds in situ, even under exuding conditions
- > Ideal structure and pore size for the recolonisation with cells and revascularisation^{6,7}
- > Full remodelling according to natural bone metabolism
- > High tolerability
- > No impairment of biomechanical properties due to thermal effects
- > Not γ -irradiated^{1,8}
- > Free of α -Gal
- > Biocompatible carrier
- > Osteoconductive (guided bone regeneration)¹⁰
- > Free of antibiotics
- > Free of preservatives
- > Can be mixed with autologous bone, blood, bone marrow and/or antibiotics⁴
- > Biocompatible carrier
- > Approved as medicinal product according to German Law on Therapeutic Products (§ 21)
- > Presentation form: putty, in glass syringe
- > Long storability of 5 years at room temperature

EXCELLENT INTRAOPERATIVE HANDLING:

- > Ready-to-use and can therefore be applied directly
- > No rehydration required
- > No thawing required
- > No mixing/stirring required
- > Does not stick to gloves

Advantage of allogeneic (human) materials

- > Allogeneic DBM demonstrates better properties than xenogeneic and synthetic materials in cell colonisation studies^{6,7}
- > Comparably high fusion rates as for the use of an autologous bone graft²

Fractures in the region of the spinal column **93.8 %** fusion rate³




Fig 2:
Example of a posterior lateral fusion with DBM pastös in an L2 fracture (10 months after stabilisation)³

DBM pastös is suitable as bone graft for PLF in traumatic thoracolumbar fractures. The use of DBM pastös in posterior lumbar fusion (PLF) prevents site morbidity. Simultaneously, the time for surgery is reduced by approximately 19 minutes¹¹. This makes DBM pastös largely cost-neutral.*

*According to a paper by Fleischer et al. in the Deutsches Ärzteblatt in 2012, every minute of surgery costs an average of EUR 45. At an inflation rate of 2.5 per cent per annum, we can therefore assume a minimum cost of EUR 50 per minute for 2017. If one multiplies this value with the minutes of surgery saved, this shows DBM pastös to be absolutely cost-neutral.


Degradation and conversion of DBM during the course of healing




The times given apply to the sheep model

- 1 0 weeks** DBM particles inserted into the defect are occupied by osteoclasts. The particles are degraded and utilised for new bone formation.
- 2 3 weeks** New bone is created at the defect margins, the first DBM particles are now occupied by osteoblasts. The particles are integrated into the new bone network.
- 3 6 weeks** During the further course, parallel bone formation and degradation (remodelling).
- 4 9 weeks** The defect has healed. The new bone matrix includes partly well integrated DBM particles.

HIGHER CELL COLONISATION WITH HUMAN DBM



Bone marrow cells adhere significantly better to human DBM than to bovine cancellous bone.⁶



After 21 days of culture, the cell count and metabolic activity on colonised DBM were significantly higher than on bovine cancellous bone and β-TCP.⁶

A safe and biocompatible carrier

Carrier properties

- > Hyaluronic acid is a polysaccharide which is formed in the human body from plasma membrane proteins.
- > Hyaluronic acid plays an important role in cell proliferation, migration, adhesion and angiogenesis.
- > Sodium hyaluronate improves the storage stability of the tissue graft.
- > DBM pastös exclusively contains high quality sodium hyaluronate pharmacopoeial quality, manufactured by fermentation in compliance with the GMP Guideline, Annex 1
- > The sodium hyaluronate used for DIZG DBM pastös is not of animal raw materials.

DBM pastös pack sizes

Article number	Grain size	Pack size
OL3550	0.212 - 0.85 mm	0.5 cm ³
OL3551	0.212 - 0.85 mm	1.0 cm ³
OL3552	0.212 - 0.85 mm	2.5 cm ³
OL3553	0.212 - 0.85 mm	5.0 cm ³
OL3554	0.212 - 0.85 mm	10.0 cm ³

REFERENCES

1. **Al Kajal, T. et al.** Evaluation of the effect of a gamma irradiated DBM-pluronic F127 composite on bone regeneration in Wistar rat. *PLoS ONE* 10, e0125110 (2015)
2. **Balling, H. et al.** Ventral, monosegmental spondylodesis of the thoracolumbar spine with demineralised bone matrix (DBM pastös) is a successful alternative to autogenous iliac crest bone graft – a prospective randomised clinical trial. *German Congress for Orthopaedics and Accident Surgery* (2012)
3. **Baumann, F. et al.** Posterolateral fusion in acute traumatic thoracolumbar fractures: a comparison of demineralized bone matrix and autologous bone graft. *Acta Ghir Orthop Traumatol Cech* 82, 119–125 (2015)
4. **Bormann, N. et al.** Analysis of parameters influencing the release of antibiotics mixed with bone grafting material using a reliable mixing procedure. *Bone* 59, 162–172 (2014)
5. **Ghanaati, S. et al.** Potential lack of 'standardized' processing techniques for production of allogeneic and xenogeneic bone blocks for application in humans. *Acta Biomater* 10, 3557–3562 (2014)
6. **Henrich, D. et al.** Characterization of bone marrow mononuclear cells on biomaterials for bone tissue engineering in vitro. *BioMed Research International* 2015, 762407–12 (2015)
7. **Janko, M. et al.** Comparison of three different types of scaffolds preseeded with human bone marrow mononuclear cells on the bone healing in a femoral critical size defect model of the athymic rat. *J Tissue Eng Regen Med* 1–37 doi:10.1002/term.2484 (2017)
8. **Kennedy, J. F. et al.** Sterilisation of tissues using ionizing radiations. *CRC Press* (2005)
9. **Kuhls, R. et al.** Human demineralised bone matrix as a bone substitute for reconstruction of cystic defects of the lower jaw. *Cell Tissue Banking* 2, 143–153 (2001)
10. **Wildemann, B. et al.** Quantification of growth factors in allogenic bone grafts extracted with three different methods. *Cell Tissue Banking* 8, 107–114 (2006)
11. **Zundler, M.** Demineralised bone matrix as alternative to autologous cancellous bone in the dorsal fusion of traumatic thoracolumbar fractures and pseudo-arthroses. *Dissertation, Med. Faculty University Clinic Regensburg* 1–125 (2015)

Name of the medicinal product:	Demineralised human bone matrix, lyophilised, DIZG
Marketing authorisation number:	PEI.H.03358.01.1
Composition:	Parts of decalcified human bone tissues. DBM pastös contains: 30% demineralised bone matrix (DBM) in a mixture of sodium hyaluronate, sodium chloride, sodium monophosphate and sodium diphosphate in water for injection.
Pharmaceutical forms:	in putty form
Therapeutic indications:	for implantation. For filling bone defects in various surgical specialist disciplines.
Contraindications:	Use in necrotic host sites is contraindicated. Should be used only when strictly indicated in poorly perfused or infected host sites, due to the poorer healing rate.
Undesirable effects:	None known
Sales limitations:	prescription only

SCANDINAVIAN DISTRIBUTOR

PUREMED APS

Universitetsparken 7
4000 Roskilde
Danmark

Phone +45 31 31 19 25
info@puremed.dk

www.puremed.dk

DIZG Deinkl. momsutsches Institut für Zell- und Gewebeersatz Gemeinnützige Gesellschaft mbH

Innovationspark Wuhlheide
Köpenicker Straße 325
D-12555 Berlin

Phone +49 30 5770 7806 0
Fax +49 30 6576 3055
distribution@dizg.de

www.dizg.de



DBMB01/2018EN